

Thin Films

The preferences of our thin films:

- 25% more pre-cut films per package than other brands
- dust proof anti static packaging
- lot number tracking system
- rolls, pre-cut circles & sheets



Item No.	Description	Thickness μm	Packaging
TF-025	Microporous		Roll 6,4cm x 5,1m
TF-050	Microporous Teflon		Roll 6,4cm x 5,1m
TF-112	Mylar®	12	Roll 7,6cm x 91,4m
TF-115	Mylar®	1,5	Roll 7,6cm x 91,4m
TF-115-345	Mylar®	1,5	Square 7,6cm x 7,6cm (500 pcs.)
TF-125	Mylar®	2,5	Roll 7,6cm x 91,4m
TF-125-255	Mylar®	2,5	Circle 6,4cm (500 pcs.)
TF-125-345	Mylar®	2,5	Square 7,6cm x 7,6cm (500 pcs.)
TF-12560	PolyRep® (replacement for polycarbonate)	Mylar® 2,5 + PP 6,0	2 Rolls a 7,6cm x 91,4m
TF-135	Mylar®	3,5	Roll 7,6cm x 91,4m
TF-135-255	Mylar®	3,5	Circle 6,4cm (500 pcs.)
TF-135-355	Mylar®	3,5	Circle 8,9cm (500 pcs.)
TF-135-345	Mylar®	3,5	Square 7,6cm x 7,6cm (500 pcs.)
TF-160	Mylar®	6,0	Roll 7,6cm x 91,4m
TF-160-255	Mylar®	6,0	Circle 6,4cm (500 pcs.)
TF-160-355	Mylar®	6,0	Circle 8,9cm (500 pcs.)
TF-160-345	Mylar®	6,0	Square 7,6cm x 7,6cm (500 pcs.)
TF-160-4	Mylar®	6,0	Square 10,16cm x 91,4m
TF-212	Polypropylene	12,0	Roll 7,6cm x 91,4m
TF-212-345	Polypropylene	12,0	Square 7,6cm x 7,6cm (500 pcs.)
TF-240	Polypropylene	4,0	Roll 7,6cm x 91,4m
TF-240-255	Polypropylene	4,0	Circle 6,4cm (500 pcs.)
TF-240-30	Polypropylene	4,0	Circle 7,6cm (100 pcs.)
TF-240-60	Polypropylene	4,0	Circle 6,4cm (60 pcs.)
TF-240-345	Polypropylene	4,0	Square 7,6cm x 7,6cm (500 pcs.)
TF-240-355	Polypropylene	4,0	Circle 8,9cm (500 pcs.)
TF-250	Polypropylene	5,0	Roll 7,6cm x 91,4m
TF-250-345	Polypropylene	5,0	Square 7,6cm x 7,6cm (500 pcs.)
TF-260	Polypropylene	6,0	Roll 7,6cm x 91,4m
TF-260-255	Polypropylene	6,0	Circle 6,4cm (500 pcs.)
TF-260-355	Polypropylene	6,0	Circle 8,9cm (500 pcs.)
TF-260-345	Polypropylene	6,0	Square 7,6cm x 7,6cm (500 pcs.)
TF-475	Kapton®	7,5	Roll 7,6cm x 15,2m
TF-475-345	Kapton	7,5	Square 7,6cm x 7,6cm (500 pcs.)
TF-500	XRF Tape	-	Roll 25mm x 66m
TF-LS-240-255	Polypropylene Low Sulfur	4,0	Circle 6,4cm (500 pcs.)

