

RackLock DigiTUBEs

Introducing the latest design for disposable digestion tubes

- * RackLock feature allows one-handed closure of screw cap digestion tubes
- * Permits dilution to volume and tube closure in an easy two-handed operation
- * Saves sample preparation time
- * Accurate 50 ml & 100ml graduations meeting Class A specifications for dilutions
- * Ultra low leachable metal content to eliminate contamination
- * Leak-free screw cap to prevent spills
- * Used in *DigiPREP* or other graphite block systems



The revolutionary *DigiTUBE* replaces:



- * Volumetric flasks
- * Graduated cylinders
- * Autosampler tubes
- * Storage containers
- * Digestion tubes

This single use, all-in-one universal vessel saves time, money, and eliminates the need to stock 5 different vessels!

Also Available for the Cleanroom

DigiTUBEs are available in shrinkwrap packaging to eliminate contamination from storage and handling.

USA

Tel.: (800) 361-6820
Fax: (800) 253-5549

Canada / International

Tel.: (800) 361-6820 /+1 (514) 457-0701
Fax: (800) 253-5549 /+1 (514) 457-4499

Europe

Tel.: +33 (0) 1 69 18 71 17
Fax: +33 (0) 1 60 92 05 67

China

Tel.: +86 (10) 58032301
Fax: +86 (10) 58032302

50ml DigiTUBES & Accessories

Catalogue No.	Description	Quantity
010-500-261	50ml <i>DigiTUBE</i> with RackLock, blue caps	750
010-500-262	50ml <i>DigiTUBE</i> with RackLock, no caps	750
010-500-263	50ml <i>DigiTUBE</i> without RackLock, blue caps	750
010-500-264	50ml <i>DigiTUBE</i> without RackLock, no caps	750
010-500-362	50ml <i>DigiTUBE</i> with RackLock, blue caps, shrink wrapped, <i>HandyPAK</i>	30 x 25
010-500-364	50ml <i>DigiTUBE</i> without RackLock, blue caps, shrink wrapped, <i>HandyPAK</i>	30 x 25
010-500-461	50ml <i>DigiTUBE</i> with RackLock, blue caps, foam rack, shrink wrapped - <i>HandyPAK Plus</i>	24 x 24
010-500-463	50ml <i>DigiTUBE</i> without RackLock, blue caps, foam rack, shrink wrapped - <i>HandyPAK Plus</i>	24 x 24
010-500-021	50ml <i>DigiTUBE</i> 24 Position RackLock Rack, for use with <i>DigiPREP</i> (48 positions)	1
010-500-025	50ml <i>DigiTUBE</i> 18 Position RackLock Rack, for use with HotBlock™	1
010-500-026	50ml <i>DigiTUBE</i> Plastic Storage Racks	5
010-500-081	50ml <i>DigiTUBE</i> Disposable Polypropylene Watch Glasses	1000
010-505-021	50ml <i>DigiTUBE</i> 24 position RackLock Digestion Rack, for use with <i>DigiPREP Jr.</i> (24 positions)	1

Order *DigiTUBES* without caps and customize 3 colors per box!

010-500-060	Blue Caps for <i>DigiTUBES</i>	250
010-500-130	Yellow Caps for <i>DigiTUBES</i>	250
010-500-140	Orange Caps for <i>DigiTUBES</i>	250
010-500-150	Red Caps for <i>DigiTUBES</i>	250
010-500-172	Clear Caps for <i>DigiTUBES</i>	250

100ml DigiTUBES & Accessories

Catalogue No.	Description	Quantity
010-501-261	100ml <i>DigiTUBE</i> with RackLock, blue caps	300
010-501-262	100ml <i>DigiTUBE</i> with RackLock, no caps	300
010-501-263	100ml <i>DigiTUBE</i> without RackLock, blue caps	300
010-501-264	100ml <i>DigiTUBE</i> without RackLock, no caps	300
010-501-028	100ml <i>DigiTUBE</i> Foam Storage Racks	5
010-501-060	Blue caps for <i>DigiTUBES</i>	100
010-501-081	100ml <i>DigiTUBE</i> Disposable Polypropylene Watch Glasses	500
010-501-510	Digestion Rack for <i>DigiPREP JR.</i> , 100ml	1
010-501-530	Digestion Rack for <i>DigiPREP MS.</i> , 100ml	1
010-501-550	Digestion Rack, for <i>DigiPREP LS.</i> , 100ml	1



Did you know about the *DigiTUBE* Calibration Certificate?

Every case of *DigiTUBES* comes complete with a lot specific Calibration Certificate that guarantees Class A tolerance at the 50 ml and 100 ml* line of every *DigiTUBE*. Quality, consistency, and cost effectiveness... that's the **SCP SCIENCE** way.

*for 100ml *DigiTUBE*

USA

Tel.: (800) 361-6820
Fax: (800) 253-5549

Canada / International

Tel.: (800) 361-6820 /+1 (514) 457-0701
Fax: (800) 253-5549 /+1 (514) 457-4499

Europe

Tel.: +33 (0) 1 69 18 71 17
Fax: +33 (0) 1 60 92 05 67

China

Tel.: +86 (10) 58032301
Fax: +86 (10) 58032302