

COMECTA Infrared thermometer "HT-100"

OPTIC
ivymen
SYSTEM

NEW

Non-contact infrared thermometer to take body temperature.
Measurement option in dual mode "Front / object", measurement in 1 second.
F1: Temperature alarm point.
F2: Temperature compensation.
F3: Open / closed warning tone.
F4: Body / object temperature.

Part No. **5903028**

Measurement distance: 5-20 cm.

Operating environment: 10-40°C ≤80% humidity.

Object: 0-100 °C.

Body mode: 32-42.9 °C.

Memory: 32 sets.

Accuracy: ± 0.2 °C.

Screen resolution: 0.1 °C.

Automatic shutdown when idle.

Unit size: 88x44x61 mm.

Weight: 160g



QUANTANALITICA

Via Guido Rossa 14
23875 Osnago (LC) - Italy
Tel. +39 039 9280061
Fax. +39 039 9289636
info@quantanalitica.com
www.quantanalitica.com