

Lithium Borates

Borate flux with guaranteed purity, **especially designed for use with the BORAMAT®**. The flux consists of micro beads allowing precise dosing and creating a perfect surface when dissolved in a fusion machine. It is anhydrous and has a low loss on ignition.

The BORAMAT® has to be operated with the micro bead flux listed below.



| Item No. | Description | Quantity in kg | Packaging |
|-------------|---|----------------|-----------|
| FX-X65-2P | Lithium tetraborate: Lithium metaborate 66:34 | 2 | Bottle |
| FX-X65-25P | Lithium tetraborate: Lithium metaborate 66:34 | 25 | Barrel |
| FX-X100-2P | Lithium tetraborate | 2 | Bottle |
| FX-X100-25P | Lithium tetraborate | 25 | Barrel |

Additional mixtures and additives on request.

