

FLUXANA[®]

XRF Application Solutions

Liquid Analysis Technology

Films

- Mylar, Polypropylene, Kapton

Reference Materials

Preparation Tool

Grease Analysis



Thin Films

The **advantages** of our thin films:

- 25% more pre-cut films per package than other brands
- dust proof, antistatic packaging
- lot number tracking system
- rolls, pre-cut circles & sheets



Item Nr.	Description	Thickness μm	Packaging
TF-025	Microporous		Roll 6.4 cm x 5.0 m
TF-050	Microporous Teflon		Roll 6.4 cm x 5.0 m
TF-112	Mylar®	12	Roll 7.6 cm x 91.4 m
TF-115	Mylar®	1.5	Roll 7.6 cm x 91.4 m
TF-115-345	Mylar®	1.5	Square 7.6 cm x 7.6 cm (500 pcs.)
TF-125	Mylar®	2.5	Roll 7.6 cm x 91.4 m
TF-125-255	Mylar®	2.5	Circle 6.4 cm (500 pcs.)
TF-125-345	Mylar®	2.5	Square 7.6 cm x 7.6 cm (500 pcs.)
TF-12560	PolyRep® (replacement for polycarbonate)	Mylar® 2.5 + PP 6.0	2 Rolls a 7.6 cm x 91.4 m
TF-135	Mylar®	3.5	Roll 7.6 cm x 91.4 m
TF-135-255	Mylar®	3.5	Circle 6.4 cm (500 pcs.)
TF-135-355	Mylar®	3.5	Circle 8.9 cm (500 pcs.)
TF-135-345	Mylar®	3.5	Square 7.6 cm x 7.6 cm (500 pcs.)
TF-160	Mylar®	6.0	Roll 7.6 cm x 91.4 m
TF-160-255	Mylar®	6.0	Circle 6.4 cm (500 pcs.)
TF-160-355	Mylar®	6.0	Circle 8.9 cm (500 pcs.)
TF-160-345	Mylar®	6.0	Square 7.6 cm x 7.6 cm (500 pcs.)
TF-160-4	Mylar®	6.0	Roll 10.1 cm x 91.4 m
TF-212	Polypropylene	12.0	Roll 7.6 cm x 91.4 m
TF-212-345	Polypropylene	12.0	Square 7.6 cm x 7.6 cm (500 pcs.)
TF-240	Polypropylene	4.0	Roll 7.6 cm x 91.4 m
TF-240-255	Polypropylene	4.0	Circle 6.4 cm (500 pcs.)
TF-240-30	Polypropylene	4.0	Circle 7.6 cm (100 pcs.)
TF-240-60	Polypropylene	4.0	Circle 6.4 cm (60 pcs.)
TF-240-345	Polypropylene	4.0	Square 7.6 cm x 7.6 cm (500 pcs.)
TF-240-355	Polypropylene	4.0	Circle 8.9 cm (500 pcs.)
TF-250	Polypropylene	5.0	Roll 7.6 cm x 91.4 m
TF-250-345	Polypropylene	5.0	Square 7.6 cm x 7.6 cm (500 pcs.)
TF-260	Polypropylene	6.0	Roll 7.6 cm x 91.4 m
TF-260-255	Polypropylene	6.0	Circle 6.4 cm (500 pcs.)
TF-260-355	Polypropylene	6.0	Circle 8.9 cm (500 pcs.)
TF-260-345	Polypropylene	6.0	Square 7.6 cm x 7.6 cm (500 pcs.)
TF-412	Kapton®	12.7	Roll 7.6 cm x 15.2 m
TF-475	Kapton®	7.5	Roll 7.6 cm x 15.2 m
TF-475-345	Kapton®	7.5	Square 7.6 cm x 7.6 cm (500 pcs.)
TF-500	XRF Tape	-	Roll 25 mm x 66 m
TF-LS-240-255	Polypropylene Low Sulfur	4.0	Circle 6.4 cm (500 pcs.)



FLUXANA® GmbH & Co. KG
Borschelstraße 3, 47551 Bedburg-Hau, Germany
Tel.: +49 (0) 2821 997 32-0
Fax: +49 (0) 2821 997 32 29
E-Mail: info@fluxana.de
Web: www.fluxana.com

Amtsgericht Kieve: HR-A 2935, HR-B 8211
Ust-IdNr.: DE 814692564, Steuer-Nr. 116/8755/0442
Finanzamt Kieve


QUANTANALITICA

Via Guido Rossa 14
23875 Osnago (LC) - Italy
Tel. +39 039 9280061
Fax. +39 039 9289636
info@quantanalitica.com
www.quantanalitica.com

Official agent